



CUSTOM BRACING

Ossur Unloader One



- Osteoarthritis
- Meniscal Tear

Used for

Ossur Function

- Pull mechanism opens up the affected compartment, allowing for improved pain relief with less joint compression.
- Lightweight and easy to put on and adjust tension.

VQ Catalyst Propel OA



- Osteoarthritis (Propel OA)
- ACL MCL LCL sprain/tear

VQ Function

- Push mechanism opens up the affected compartment, allowing for increased space at the joint.
- Great for high activity levels.
- Very rigid for contact sports.

Don Joy Defiance III



- Osteoarthritis
- ACL MCL LCL sprain/tear

Don Joy Function

- Custom carbon fibre rigid frame to provide superior stability and support to surrounding knee ligaments.
- unloader option for OA or meniscus tears to decrease pressure on affected side.
- Strongest of all DonJoy knee braces – can be used for contact or non contact sports.

DonJoy ROAM OA



- Unicompartment OA

DonJoy ROAM OA

- Light weight low profile single hinge brace.
- Provides mobility and pain relief by unloading the affected compartment.
- Easy application and adjustable tension.

Spring loaded



- Mild Osteoarthritis
- Patellofemoral OA

Spring loadedFunction

- Spring-loaded mechanism assists with the extension of the knee.



226-336-8961
EXT 2



201- 245 The Boardwalk Waterloo
ON. N2T0A6



THEARMOURYCLINIC.CA